

2009 SPE-DEX 4790 Master Parts List

Service Parts, Spare Parts, Accessories and Consumables

Effective January 1, 2009

Item No.	Description
50-2498	4790 - SAMPLE COLLECT VALVE KIT (1)VALVE, SOLENOID, 3-WAY (2)TERM, M, MINI-FIT JR, AWG22-28
50-2535	REPL VALVE 4790 WATER/WASTE KIT (1)KIT, VALVE, SOL/H2O WASTE 4790 (2)TERM, M, MINI-FIT JR, AWG22-28
50-2536	REPL VALVE 4790 PREWET & RINSE KIT (1)Valve (Manifold Block) (2)TERM, M, MINI-FIT JR, AWG22-28
48-2604	4790 - THERMISTOR ASSEMBLY (KIT) (2)4790 - THERMISTOR ASSEMBLY
48-2605	SAMPLE COLLECT CHECK VALVE (KIT)
48-2606	DISPENSE RING AND TUBE ASSY (KIT) (1)RING, PRE-WET, 4790 (1) RING, DISPENSE, PRE-WET, 4789
48-2607	H2O Inlet Valve, 4790 (KIT) (1)BODY, FUNNEL, 4790 (2)LINER, VALVE, 4790 (1)NUT, SPANNER, 4790 (1)ASSY, STEM, VALVE, H2O INLET (1)PACKING, VALVE, H2O INLET,4790 (1)STEM, VALVE, H2O INLET, 4790 (2)8-32 x 5/8 PHIL MACH SCREW, SS (1)10-32 x 3/8" SCREW, SET SCKET (1)PIN, SPRING, 420 SS 5/32"x3/4" (1)SPRING, WAVE
50-0620-PL	ASSY, VALVE, H2O INLET, 4790 (2)LINER, VALVE, 4790 (1)WATER INLET VALVE STEM ASSY (1)PACKING, VALVE, H2O INLET,4790 (1)STEM, VALVE, H2O INLET, 4790
50-2457	ASSY, HOLDER, BOTTLE, 4790

(1)HOLDER, BOTTLE, 4790
(1)TUBE, RINSE, 4790
(1)O-RING, HOLDER, BOTTLE, SAMPLE
(1)Warning Label

28-0615 **MOTOR, SPRING**

99-0766 **TOOL - MOLEX PIN EXTRACTOR**

02-0727 **4790 - CHECK VALVE TOOL**

165-0712 **HYDROPHILIC SOLUTION, 250ML**
(1)Hydrophilic solution 250 ml
(3) 20ml polyethylene scintillation bottle

16-0626 **SENSOR, PHOTOINTERRUPTER, IR**

Accessories

30-0628 **O-RING, HOLDER, BOTTLE, SAMPLE**

02-0602 ADAPTOR, COLLECTION, 4790

22-0007 Nut, Compression, Peek, 5/16"

22-0008 Ferrule, Compression, 3/16", B

99-0664 C-Ring Pliers

02-0682 PLUG, LUER, F, PVDF, MODIFIED

22-0011-02 Tee, Luer, Kynar

22-0012-02 Cap, Luer Lock, Rotating, Kyna

22-0021-02 Elbow, Kynar, 1/8 x Luer

22-0047 Reducer, Brass 1/8" x 1/4" NPT

22-0049-02 CAP, PLUG, KYNAR

22-0197-02 Male Luer Lock to 1/8 Barb - K

120-4001-EX-END 4790 - SOL. MANIFOLD END BLOCK

120-4001-EX-THRU 4790 - SOL.MANIFOLD THRU BLOCK

02-0005 Solvent Delivery Bottle Plug

02-0006 Solvent Bottle Gasket 1.74x1.1

02-0562 WASTE CARBOY - STOPPER, 20L

02-0591-01 RING, PRE-WET, 4790

02-0594 BODY, FUNNEL, 4790

02-0595 LINER, VALVE, 4790

02-0597-01 HOLDER, BOTTLE, 4790

02-0601 NUT, SPANNER, 4790

02-0608 TUBE, RINSE, 4790

02-0618 CURTAIN , 4790

02-0621-01 RING, DISPENSE, PRE-WET, 4790

02-0623-03 RISER, EXTENDED, 47MM 3K/4790

02-0738 STPK Plug, 2mm bore, Teflon, W

02-0739 Support Screen, LLE - DD

02-0882	4790 - EXHAUST HOSE 10 FEET
02-0882-01	4790 - EXHAUST HOSE 1-PLACE
02-0882-02	4790 - EXHAUST HOSE 2-PLACE
02-0882-03	4790 - EXHAUST HOSE 3-PLACE
02-0882-04	4790 - EXHAUST HOSE 4-PLACE
02-0882-05	4790 - EXHAUST HOSE 5-PLACE
02-0882-06	4790 - EXHAUST HOSE 6-PLACE
02-0882-07	4790 - EXHAUST HOSE 7-PLACE
02-0882-08	4790 - EXHAUST HOSE 8-PLACE
02-0882-ADD-1	4790 - EXHAUST HOSE ADD-ON
02-1573-01	CARTRIDGE ADAPTOR - EXTENDED
50-2532	LOW VACUUM CARTRIDGE REGULATOR
05-0365	Tubing, Microline, #1861 1/8"ODx
06-0110-01	VALVE, CHECK, SAMPLE COLLECT
06-0709	VACUUM - BLEED VALUE, 1/4"
07-2399	JT BAKER SPEEDISK RESERVOIR
100-10119-CE	4790 - COMMUNICATION CABLE 7'
160-0001-02	VOA VIAL ADAPTER - 19/22, 40ML
160-0008	VOA VIAL 40ML - 72CT.
165-00JT-01	SPE Disk Holder, JT Baker , 47 (pkg of 2)
19-0654	4790 Packing Matl
20-0016	FEET, RUBBER, LARGE
22-0176	Tefzel Plug for 1/4 - 28 Fitti
22-0311	NUT, HEADLESS, 1/4-28, PPS
22-0312	FERRULE W/ RING, 1/8", TEFZEL
22-0687	Clamp, #19, Blue
27-0834-01	SEP. FUNNEL - 125ML, 19/22
40-012-HT	Inlet Solvent Filter, A-310
50-001-HT	O&G - TRANSFER PIPETTES, 250#
50-002-HT	O&G - PANS, 70MM, 55ML, 100#
50-002-02-HT	O&G - PANS, 105MM, 125ML, 100#
50-0027-03	4790 - PRESSURE REG. BRACKET
50-0043-02	4790 - SOL. VAC. MANIF. ASSY.
50-0051-02	4790 - H2O VAC. REGULATOR. KIT
50-0103-01	4790 - SOL.DELIV LINES, 3/16"
50-0153	ASSEMBLY, THERMISTOR, 4790
50-0172-01	Solvent Bottle Plug Assembly,
50-0173-04	4790 - SOL.DELIV. LINES, 3/32"
50-0211-02	4790, 1/3000XL - VAC. LINE KIT
50-0225	4790 - SOL. & H2O WASTE LINES
50-0632-02	ASSEMBLY, WHEEL, PIVOT, 4790
50-1686-00	ASSY 20L WW CARBOY W/ DRAIN
50-2111	4790 - DUAL PH KIT, 38X430 BOT
EZ-SEALS-HT	4790-AI SQUARES, W/CASE, 500#
27-0042	BOTTLE - PLAS. COAT., 2.5L
27-0271	BOTTLE - BOS. ROUND - CLEAR 1L
27-0271-01	BOTTLE - BOS. ROUND - 1L (Case of 12)

27-0476 ERLENMEYER - 125ML, 19/22

Cap Adapters

150-0001 33 x 400 Adapter, 1L, Boston Round
150-0008-01 33 x 430 Adapter, 1L Narrow Mouth
150-0023 38X400 Adapter
150-0021 38X430 Adapter
150-0019 43 x 400 Adapter, 250 ml bottle
150-0004-01 43MM, NALG.CAP ADAP.
150-0020 45 x 400 Adapter, 250 ml bottle
150-0012-02 48 x 400 Adapter
150-0009-05 53 x 400 Adapter
150-0010 58 x 400 Adapter
150-0014 63 x 400 Adapter
150-0005 70 x 400 Adapter, 1L bottle
150-0005-02 70 x 400 Adapter, Ball/Kerr Jar
150-0015 83 x 400 Adapter
150-0015-02 83X400, METAL Adapter
150-0006 89 x 400 Adapter
150-2260 CAP ADPTR ASSY 45MM SCHOTT
150-2424 WINCHESTER BOTTLE CAP ADAPTOR
150-0024 4790 - CAP ADAPTER 40ML VOA
27-0003-03 CAP ADAPTER - CAP ONLY, 38X430

Gaskets

02-0560-02 33 mm gasket (1.255" OD x 0.97" ID)
02-0560-05 48 mm gasket (1.855" OD x 1.52" ID)
02-0560-06 53 mm gasket (2.075" OD x 1.810" ID)
02-0560-07 58 mm gasket (2.185" OD x 1.88" ID)
02-0560-08 63 mm gasket (2.445" OD x 2.06" ID)
02-0560-09 70 mm gasket (2.74" OD x 2.38" ID)
02-0560-10 83 mm gasket (3.27" OD x 2.85" ID)
02-0560-11 89 mm gasket (3.52" OD x 3.10" ID)
02-0560-12 43 mm gasket (1.72" OD x 1.118" ID)
02-0560-21 Ball Jar Gasket (3.41" OD x 2.88" ID)

Disk Holders and Accessories

50-0807-01 SPE Disk Holder, 47 mm , 4790
50-0809-01 SPE Disk Holder, 90 mm , 4790
02-0592 Disk Holder Luer Fitting - MODIFIED
110-0017 PURGE OR 47 mm, 3M Disk Oil & Grease Gasket
'02-2492 90 mm, 3M Disk Oil & Grease Gasket
31-0454 47 mm Fine Mesh Screen - (for dirty samples)
31-0457 90 mm Fine Mesh Screen - (for dirty samples)
31-0661 DISK HOLDER - 47MM, SUP.SCREEN
31-0662 DISK HOLDER - 90MM, SUP.SCREEN

Vacuum Pump

50-2480-01	GAST 120V 60HZ VAC PMP PLUMBED
50-2480-02	GAST 220V 50/60HZ VAC PMP PLMB
50-2483-01	KNF EPM 120-220V 60HZ PLUMBED
50-2483-02	KNF EPM 220V 50HZ 1PH PLUMBED
50-2483-03	KNF EPM 220V 50HZ 3PH PLUMBED

Pump Accessories

49-2481-01	GAST NEPM VAC PMP PLUMBING KIT
49-2481-02	KNF EPM VAC PUMP PLUMBING KIT

SPE Consumables

47-2346-02	ATL DISK C18 TYP B 50MM - BX24
47-2346-06	ATL DISK DVB 50MM - BX24
47-2346-05	OASIS DISK HT 1G 30U - PKG20
50-003-HT	O&G - STANDARD - 26ML, 4#
1664-47-HT	PACIFIC DISK - 1664A - 47MM
1664-90-HT	PACIFIC DISK - 1664A - 90MM
FFP-47-HT	FAST FLOW PREFILTER 47 BX 50
FFP-90-HT	FAST FLOW PREFILTER 90 BX 50
47-001-HT	3M - DISK C8, 47MM, 20#
47-002-HT	3M Empore Disks C18 47mm (Octa
47-003-HT	3M - DISK SDB-RPS, 47MM, 20#
47-004-HT	3M - DISK ANION-SR, 47MM, 20#
47-005-HT	3M - DISK O&G, 47MM, 20#
47-007-HT	3M - DISK SDB-XC, 47MM, 20#
47-008-HT	3M - DISK O&G, 90MM, 10#
47-009-HT	3M Empore 90mm Disks,C18 (Octa
47-011-HT	3M - DISK SDB-XC, 90MM, 10#
50-019-HT	JT BAKER SPD.DSK.- O&G, 20#
50-643-HT	JT BAKER SPD.DSK.- C18, 20#
50-644-HT	JT BAKER SPD.DSK.- DVB, 20#
50-645-HT	JT BAKER SPD.DSK. - SAX, 20#
50-646-HT	JT BAKER SPD.DSK.- C8, 20#
50-647-HT	JT BAKER SPD.DSK.- C18XF, 20#
50-742-HT	JTBKR.SD-DVB H2O PHI, 5G, 20#
50-799-HT	JTBKR.SD-DVB H2O PHI, 1G, 20#
50-845-HT	JTBKR.SD-DVB H2O PHI, 8270,20#
40-001-HT	3M - FILTERAID, 1.5KG - FA400